

Product Data Sheet

Product Name: E3 ligase Ligand-Linker Conjugates 33 Hydrochloride

Cat. No.: GC35938

Chemical Properties

Cas. No. 1239866-59-1

SMILES NCCOCCOC([C@@H](NC([C@H](C(CC1=CC=CC=C1)NC(OCC2C3=C(C=CC=C3)C4=C2C=CC=C4)=O)O)=O)CC(C)C)=O.Cl


Formula $C_{35}H_{44}ClN_3O_7$ M.Wt 654.19

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

E3 ligase Ligand-Linker Conjugates 33 Hydrochloride incorporates a cIAP ligand for the E3 ubiquitin ligase, and a PROTAC linker. E3 ligase Ligand-Linker Conjugates 33 can be used to design

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

PROTAC degrader[1]. cIAP1

[1]. Itoh Y, et al. Design, synthesis and biological evaluation of nuclear receptor-degradation inducers. Bioorg Med Chem. 2011 Nov 15;19(22):6768-78.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA