

Product Data Sheet

Product Name: Manganese(II) carbonate hydrate, 44-46% Mn basis (KT)

Cat. No.: GB59260

Chemical Properties

Cas. No. 34156-69-9

Formula $\text{MnCO}_3 \cdot x\text{H}_2\text{O}$

M.Wt 114.95 (anhydrous basis)


Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA