

Product Data Sheet

Product Name: (2R,2R)-PF-07258669

Cat. No.: GC73162

Chemical Properties

Cas. No. 2755890-68-5


Formula $C_{25}H_{27}FN_6O_2$ M.Wt 462.52

Solubility Storage -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

(2R,2R)-PF-07258669 is a MC4R (melanocortin 4 receptor) antagonist. (2R,2R)-PF-07258669 can be used in the study of regulating appetite and energy expenditure.

References:

[1]. Cao Q, et al. Modular, automated synthesis of spirocyclic tetradronaphthridines from primary alkylamines[J]. Communications Chemistry, 2023, 6(1): 215.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA