

Product Data Sheet

Product Name: Lei-Dab 7

Cat. No.: GC50207

Chemical Properties

Cas. No. 1061556-49-7

Formula $C_{141}H_{236}N_{46}O_{39}S_6$

M.Wt 3392.06


Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

High affinity, selective KCa2.2 (SK2) channel blocker ($K_d = 3.8$ nM). Exhibits >200-fold selectivity for KCa2.2 over KCa2.1, KCa2.3, KCa3.1, IK, Kv and Kir2.1. Increases theta-burst responses and increases LTP in rat hippocampal slices in vitro. Convulsive in vivo.

Shakkottai et al (2001) Design and characterization of a highly selective peptide inhibitor of the small conductance calcium-activated K⁺ channel, SkCa2. J.Biol.Chem. 276 43145 PMID:11527975 |Kramar et al (2004) A novel mechanism for the facilitation of theta-induced long-term potentiation by brain-derived neurotrophic factor. J.Neurosci. 24 5151 PMID:15175384 |Kasumu et al (2012) Selective positive modulator of calcium-activated potassium channels exerts beneficial effects in a mouse model of spinocerebellar ataxia type 2. Chem.Biol. 19 1340 PMID:23102227

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA