

Product Data Sheet

Product Name: PE 22-28

Cat. No.: GC72194

Chemical Properties

Cas. No. 1801959-12-5


Formula $C_{35}H_{55}N_{11}O_9$ M.Wt 773.88

Solubility Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

PE 22-28 is a TREK-1 inhibitor with IC₅₀ value of 0.12 nM. PE 22-28 also is a 7 amino-acid peptide that is used as a core sequence for preparing analogs by chemical modifications and also by substitution of amino-acids. PE 22-28 can be used for the research of depression.

PE 22-28 displays good specificity and affinity for TREK-1 channel with IC₅₀ value of 0.12 nM^[1].

PE 22-28 (i.p.; 3.0µg/kg) shows antidepressant property^[1].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Alaeddine Djillani, et al. Shortened Spadin Analogs Display Better TREK-1 Inhibition, In Vivo Stability and Antidepressant Activity. Front Pharmacol. 2017 Sep 12;8:643.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA