

Product Data Sheet

Product Name: HP210
Cat. No.: GC70973

Chemical Properties

Cas. No. 728039-52-9


Formula $C_{22}H_{19}N_3O_2S_2$ M.Wt 421.54

Solubility DMSO : 125 mg/mL (296.53 mM); Storage
Need ultrasonic; Hygroscopic DMSO
has a significant impact on the
solubility of product, please use
newly opened DMSO)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

HP210 is a selective glucocorticoid receptor modulator (SGRM). HP210 can inhibit the mRNA expression of IL-1 β and IL-6. HP210 has the potential to study inflammation-related diseases.

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Hu X, et al. Discovery of novel non-steroidal selective glucocorticoid receptor modulators by structure- and IGN-based virtual screening, structural optimization, and biological evaluation. Eur J Med Chem. 2022 Jul 5;237:114382.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA