

Product Data Sheet

Product Name: Ginsenoside Rg5

Cat. No.: GC31706

Chemical Properties

Cas. No. 186763-78-0

SMILES C[C@@]([C@@]12C)(CC[C@@]3([H])C4(C)C)[C@@](C[C@@H](O)[C@]1([H])[C@@H](/C(C)=C/C/C=C(C)/C)CC2)([H])[C@]3(CC[C@@H]4O[C@@](O[C@H](CO)[C@@H](O)[C@@H]5O)([H])[C@@H]5O[C@]([C@@H]([C@@H](O)[C@@H]6O)O)([H])O[C@@H]6CO)C


Formula C₄₂H₇₀O₁₂ M.Wt 767

Solubility DMF: 10 mg/ml, DMSO: 10 mg/ml Storage Store at -20°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Protocol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

**Kinase
experiment:**

HUVECs are cultured in 24-well plates overnight. The cells are changed to serum-free M199 and incubated for 1 h. The medium is removed, and cells are incubated with fresh serum-free medium containing 0.1 μ M-50 mM Ginsenoside Rg5 at 37°C for 20 min followed by the addition of 50 μ L (1 μ Ci) of [125I]IGF-1 and then further incubated for 10 min. The medium is decanted, and cell plates are washed twice with serum-free medium. Cells are lysed in 300 μ L of 0.1 N NaOH solution containing 0.1% SDS, transferred to scintillation vials, and mixed with 1 mL of Ultima Gold mixture solution. Cell-associated [125I]IGF-1 is analyzed in a scintillation counter. The nonspecific binding is determined by coincubation with unlabeled IGF-1 (50 nM)[1].

Cell experiment:

MCF-7 (HER2-/ER+) and MDA-MB-453 (HER2+/ER-) human breast cancer cell lines are maintained using RPMI 1640 medium supplemented with 10% (vol/vol) FBS plus 100 units/mL Penicillin and Streptomycin in a 5% carbon dioxide air incubator at 37°C. Cell cytotoxicity is measured by MTT assay. Cells are seeded in 96-well tissue culture plates at the density of 0.2×10^4 cells per well with 100 μ L medium, and are allowed to become attached for 24 h. One hundred microliters of the medium with different concentrations of Ginsenoside Rg5 (e.g., 0 μ M, 25 μ M, 50 μ M, and 100 μ M) are added to each well. At indicated times, 30 μ L MTT stock solution (3 mg/mL) are added to each well. After culturing the cells at 37°C for 2 h, DMSO is added to dissolve the formazan crystals. The absorbance is read at the wavelength of 540 nm with a microplate reader[3].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Animal experiment: Mice[2]Male ICR mice (6 to 8 weeks old), weighing 25-27 g, are used. After acclimation for one week, mice are randomly assigned into 4 experimental groups with 8 mice in each group: normal control, Cisplatin control, and Cisplatin+Ginsenoside Rg5 groups (10 and 20 mg/kg, respectively). Ginsenoside Rg5 is administered intragastrically at the dose of 10 and 20 mg/kg for 10 days. On the 7th day, animals in Cisplatin control and Ginsenoside Rg5-treated groups receive a single intraperitoneal injection of Cisplatin (25 mg/kg) to induce nephrotoxicity in mice. Mice are anaesthetized with pentobarbital, subsequently sacrificed at 72 h after Cisplatin injection (Day 10). Blood samples are collected and then centrifuged at 3000 rpm to separate the serum and stored at -20 °C for determining blood urea nitrogen (BUN) and creatinine (CRE) levels.

References:

- [1]. Cho YL, et al. Specific activation of insulin-like growth factor-1 receptor by ginsenoside Rg5 promotes angiogenesis and vasorelaxation. J Biol Chem. 2015 Jan 2;290(1):467-77.
- [2]. Li W, et al. Ginsenoside Rg5 Ameliorates Cisplatin-Induced Nephrotoxicity in Mice through Inhibition of Inflammation, Oxidative Stress, and Apoptosis. Nutrients. 2016 Sep 13;8(9). pii: E566.
- [3]. Kim SJ, et al. Anti-breast cancer activity of Fine Black ginseng (Panax ginseng Meyer) and ginsenoside Rg5. J Ginseng Res. 2015 Apr;39(2):125-34.

Background

Ginsenoside Rg5 is the main component of Red ginseng. Ginsenoside blocks binding of IGF-1 to its receptor with an IC₅₀ of ~90 nM. Ginsenoside Rg5 also inhibits the mRNA expression of COX-2 via suppression of the DNA binding activities of NF-κB p65.

Ginsenoside Rg5 plays a novel role as an IGF-1R agonist. Ginsenoside Rg5 binds to IGF-1R with an IC₅₀ value of ~90 and its angiogenic activity is inhibited by IGF-1R knockdown. To investigate the possible interaction of Ginsenoside Rg5 with IGF-1R, a docking analysis is performed. Docking results show that Ginsenoside Rg5 binds strongly at two sites, A and B, with K_d values of 20 and 27 nM, respectively, to the cysteine-rich domain of IGF-1R. Pretreatment with Rg5 blocks the binding of radiolabeled IGF-1 to HUVECs with an IC₅₀ value of ~90 μM, which is greater than an IC₅₀ value of

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

~1.4 nM for unlabeled IGF-1[1]. The results from MTT assay show that MCF-7 cell proliferation is inhibited by Ginsenoside Rg5 treatment for 24, 48 and 72 h in a dose-dependent manner. Ginsenoside Rg5 at different concentrations (0, 25, 50 and 100 μ M), induce cell cycle arrest in G0/G1 phase through regulation of cell cycle-related proteins in MCF-7 cells[3].

Ginsenoside Rg5 inhibits the mRNA expression of COX-2 via suppression of the DNA binding activities of NF- κ B p65 in lipopolysaccharides (LPS)-stimulated BV2 microglial cells. Rg5 pretreated group mice show declined expression of NF- κ B p65 and COX-2. In the group treated with low dose of Ginsenoside Rg5 (10 mg/kg), there is remarkable tubular damage and infiltration of inflammatory cells. However, at the higher dose of Ginsenoside Rg5 (20 mg/kg), tubules markedly appear histologically normal and no inflammation and cast formation is observed in kidney tissues[2].

[1]. Cho YL, et al. Specific activation of insulin-like growth factor-1 receptor by ginsenoside Rg5 promotes angiogenesis and vasorelaxation. J Biol Chem. 2015 Jan 2;290(1):467-77. [2]. Li W, et al. Ginsenoside Rg5 Ameliorates Cisplatin-Induced Nephrotoxicity in Mice through Inhibition of Inflammation, Oxidative Stress, and Apoptosis. Nutrients. 2016 Sep 13;8(9). pii: E566. [3]. Kim SJ, et al. Anti-breast cancer activity of Fine Black ginseng (Panax ginseng Meyer) and ginsenoside Rg5. J Ginseng Res. 2015 Apr;39(2):125-34.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA