

Product Data Sheet

Product Name: (3-(6-Phenylhexyl)phenyl)boronic acid

Cat. No.: GF52734

Chemical Properties

Cas. No. 1795440-43-5

Formula C₁₈H₂₃BO₂

M.Wt 282.19

Solubility

Storage Store at RT, away from moisture sealed storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA