

Product Data Sheet

Product Name: (Rac)-Atropine-d3

Cat. No.: GC68414

Chemical Properties

Cas. No. 1276197-36-4

Formula $C_{17}H_{20}D_3NO_3$ M.Wt 292.39

Solubility DMSO : 50 mg/mL (171.00 mM; Need ultrasonic); Ethanol : ≥ 16 mg/mL (54.72 mM); DMSO : ≥ 10 mg/mL (34.20 mM); PBS (pH 7.2) : ≥ 10 mg/mL (34.20 mM); DMF : ≥ 2 mg/mL (6.84 mM)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA