

Product Data Sheet

Product Name: MBHA Resin

Cat. No.: GA11267

Chemical Properties

Cas. No. N/A

Chemical Name phenyl(p-tolyl)methanamine compound with 1,4-divinylbenzene and styrene (1:1:1) hydr ochloride

SMILES CC1=CC=C(C(N)C2=CC=CC=C2)C=C1.C=CC3=CC=CC=C3.C=CC4=CC=C(C=C)C=C4.Cl

Formula $C_{32}H_{34}ClN$

M.Wt 468.07


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

MBHA Resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA