

Product Name: VHL-IN-1
Cat. No.: GC73611

Cas. No. 3033117-53-9


Formula	C ₂₈ H ₃₇ FN ₄ O ₄ S	M.Wt	544.68
---------	--	------	--------

Solubility	DMSO : 100 mg/mL (183.59 mM; Need ultrasonic)	Storage	-20°C, sealed storage, away from moisture and light
------------	--	---------	--

General tips	For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.
--------------	---

Shipping Condition	Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
--------------------	--

Structure



VHL-IN-1 (compound 30) is a ubiquitin E3 ligase von Hippel-Lindau (VHL) inhibitor (dissociation constant $K_d=37$ nM) that stabilizes and induces HIF-1 α transcriptional activity. VHL-IN-1 has potential as a HIF-1 α stabilizer and degrader of proteolytically targeted chimeras (PROTACs).

[1]. Vu LP, et al. Expanding the Structural Diversity at the Phenylene Core of Ligands for the von Hippel-Lindau E3 Ubiquitin Ligase: Development of Highly Potent poxia-Inducible Factor-1 α

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Stabilizers. J Med Chem. 2023 Sep 14.. .

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA