

Product Data Sheet

Product Name: ATI-2341 TFA

Cat. No.: GC38541

Chemical Properties

Cas. No.


Formula $C_{106}H_{179}F_3N_{26}O_{27}S_2$ M.Wt 2370.84

Solubility DMSO: 100 mg/mL (42.18 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

ATI-2341 is a potent and functionally selective allosteric agonist of C-X-C chemokine receptor type 4 (CXCR4), which functions as a biased ligand, favoring $G_{\alpha i}$ activation over $G_{\alpha 13}$. ATI-2341 activates the inhibitory heterotrimeric G protein (G_i) to promote inhibition of cAMP production and induce calcium mobilization. ATI-2341 is a potent and efficacious mobilizer of bone marrow polymorphonuclear neutrophils (PMNs) and hematopoietic stem and progenitor cells (HSPCs)[1][2].

CXCR4[1]

[1]. Quoyer J, et al. Pepducin targeting the C-X-C chemokine receptor type 4 acts as a biased agonist

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

favoring activation of the inhibitory G protein. Proc Natl Acad Sci U S A. 2013 Dec 24;110(52):E5088-97. [2]. Tchernychev B, et al. Discovery of a CXCR4 agonist pepducin that mobilizes bone marrow hematopoietic cells. Proc Natl Acad Sci U S A. 2010 Dec 21;107(51):22255-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA