

Product Data Sheet

Product Name: RIPK2 ligand 3

Cat. No.: GC80878

Chemical Properties

Cas. No. 2027561-86-8

Formula $C_{20}H_{19}N_3O_3S_2$

M.Wt 413.51


Solubility

Storage Store at 4°C, protect from light,
stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

RIPK2 ligand 3 is a RIPK2 ligand and can be used for the synthesis of PROTACs, such as PROTAC RIPK degrader-2 [1].

References:

[1]. Bondeson DP, et al. Catalytic in vivo protein knockdown by small-molecule PROTACs. Nat Chem Biol. 2015 Aug;11(8):611-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA