

Product Data Sheet

Product Name: GLP-1 receptor agonist 2

Cat. No.: GC36147

Chemical Properties

Cas. No. 2230197-64-3

SMILES O=C(C1=CC=C2N=C(CN3[C@@H](C)CN(C4=NC(OCC5=CC=C(Cl)C=C5F)=CC=C4)CC3)N(C[C@H]6OCC6)C2=C1)O


Formula $C_{30}H_{31}ClFN_5O_4$ M.Wt 580.05

Solubility DMSO: ≥ 125 mg/mL (215.50 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

GLP-1 receptor agonist 2 is a glucagon-like peptide-1 receptor (GLP-1R) agonist. GLP-1R[1].

[1]. WO 2018/109607 A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA