

Product Data Sheet

Product Name: MS102
Cat. No.: GC81468

Chemical Properties

Cas. No. 3058587-77-9


Formula $C_{24}H_{17}BrF_3N_3O_3S_2$ M.Wt 596.44

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

MS102 is an orally active ubiquitin specific peptidase 2 (USP2) inhibitor with an IC_{50} of 5.46 μM . MS102 has viable antiviral activity against ACE2 -dependent coronaviruses. MS102 significantly reduces V-domain Ig suppressor of T cell activation (VISTA) protein abundance in vitro and in vivo. MS102 can be used for the study of SARS-CoV-2. MS102 can be used in combination with anti-PD-1 immunotherapy to enhance the anti-tumor immune response [1] [2].

In Vivo, MS102 (25 mg/kg, i.p., daily for 4 days, followed by 6 hours challenge with SARS-CoV-2 WT) reduces SARS-CoV-2 viral load and lung pathology in K18-hACE2 transgenic mice^[1]. MS102

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

(10 mg/kg, intraperitoneal daily, for 21 days) enhances anti-PD-1 immunotherapy by reducing VISTA and promoting T cell activity in syngeneic mouse tumor models^[2].

In Vitro, MS102 (4-10 μ M; 24-48 h) demonstrates improved efficacy in blocking the infection of primary human bronchial epithelial (NHBE) cells by the authentic SARS-CoV-2 wild-type. MS102 has minimal cytotoxicity in HepG2 and Calu-3 cells^[1]. MS102 is effective in inhibiting the ACE2-dependent coronavirus hCoV-NL63 (an ACE2-dependent coronavirus) with an EC₅₀ of 80 nM^[1]. MS102 demonstrates significant reduction in VISTA protein levels in MDA-MB-231, SUM159, MDA-MB-436, EMT6, and RENCA cells^[2].

References:

[1]. Dang F, et al. USP2 inhibition prevents infection with ACE2-dependent coronaviruses in vitro and is protective against SARS-CoV-2 in mice. *Sci Transl Med.* 2023;15(725):eadh7668.

[2]. Li Chen, et al. Targeted destruction of VISTA boosts anti-tumor immunotherapy. *Cell Res.* 2025 Dec;35(12):987-1002.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA