

Product Data Sheet

Product Name: CRBN ligand-896

Cat. No.: GC81009

Chemical Properties

Cas. No. 2504234-93-7

Formula $C_{16}H_{11}BrN_2O_3$

M.Wt 359.17


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

CRBN ligand-896 (CRBN-1) is a E3 ligase ligand of PROTAC that can be used to synthesize PROTAC degraders targeting the epigenetic coactivator p300 [1].

References:

[1]. Abbineni C, et al. Preparation of substituted heterocyclic compounds and their derivatives as CBP and/or p300 degraders. World Intellectual Property Organization, WO2025202979 A1 2025-10-02.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA