

Product Data Sheet

Product Name: Ethyleneglycolmonopentylether

Cat. No.: GA21430

Chemical Properties

Cas. No. 6196-58-3

Formula $C_7H_{16}O_2$

M.Wt 132.2

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

The tenside C5E1 solubilizes membrane proteins.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA