

Product Data Sheet

Product Name: Sodium tetraborate pentahydrate

Cat. No.: GB61147

Chemical Properties

Cas. No. 12179-04-3

Formula $B_4H_{10}Na_2O_{12}$

M.Wt 291.3

Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA