

Product Data Sheet

Product Name: p-Tolualdehyde

Cat. No.: GC39789

Chemical Properties

Cas. No. 104-87-0

SMILES CC1=CC=C(C=O)C=C1


Formula C8H8O M.Wt 120.15

Solubility DMSO : 100 mg/mL (832.29 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

p-Tolualdehyde is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA