

Product Data Sheet

Product Name: Biotinyl-Atrial Natriuretic Factor (1-28) (human)

Cat. No.: GA20828

Chemical Properties

Cas. No. 1815618-06-4

Formula $C_{137}H_{217}N_{47}O_{41}S_4$

M.Wt 3306.75

Solubility H_2O : 5 mg/mL (1.51 mM; Need ultrasonic)

Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA