

Product Data Sheet

Product Name: Propargyl-PEG1-NHS ester

Cat. No.: GC64251

Chemical Properties

Cas. No. 1174157-65-3

Formula C₁₀H₁₁NO₅

M.Wt 225.2


Solubility DMSO : 100 mg/mL (444.05 mM;
Need ultrasonic)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Propargyl-PEG1-NHS ester is a noncleavable 1-unit PEG linker for antibody-drug-conjugation (ADC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA