

Product Data Sheet

Product Name: Mopivabil

Cat. No.: GC64604

Chemical Properties

Cas. No. 2988562-84-9

Formula C₁₄H₂₀O₃

M.Wt 236.31


Solubility DMSO : 100 mg/mL (423.17 mM;
Need ultrasonic)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Background

Mopivabil is the antagonist of angiotensin II receptor[1].

[1]. WHO Drug Information, Vol. 35, No. 4, 2021. Geneva: World Health Organization; 2022.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA