

Product Data Sheet

Product Name: immunoglobulin light chain variable region fragment [Homo sapiens]/[Mus musculus]
Cat. No.: GP10133

Chemical Properties

Cas. No.

SMILES NC(CC(C)C)C(NC(CC(C)C)C(NC(C(C)CC)C(NC(CC1=CC=C(O)C=C1)C(NC(CC2=CNC3=C2C=CC=C3)C(NC(C)C(NC(CO)C(NC(C(O)C)C(NC(CCCNC(N)=N)C(O)=O)=O)=O)=O)=O)=O)=O


Formula $C_{54}H_{83}N_{13}O_{13}$ M.Wt 1122.32

Solubility $\geq 112.2\text{mg/mL}$ in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA