

Product Data Sheet

Product Name: 2,4-Dimethoxyprimidine-5-boronic acid(contains varying amounts of Anhydride)
Cat. No.: GD08905

Chemical Properties

Cas. No. 89641-18-9

Formula C₆H₉BN₂O₄ M.Wt 183.96

Solubility DMSO : 100 mg/mL (543.60 mM; Storage 2-8°C argon filling
Need ultrasonic)

General For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the
tips ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT , or
Condition blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Background

2,4-Dimethoxypyrimidine-5-boronic acid is a biochemical reagent that can be used as a biological material or organic compound for life science related research.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA