

Product Data Sheet

Product Name: Methacrylate C6 phosphoramidite

Cat. No.: GC81465

Chemical Properties

Cas. No. 1226983-36-3


Formula $C_{19}H_{36}N_3O_3P$ M.Wt 385.49

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Methacrylate C6 phosphoramidite is a phosphoramidite that can be used in the synthesis of oligonucleotides.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA