

Product Data Sheet

Product Name: PROTAC STAT6 degrader-3

Cat. No.: GC80949

Chemical Properties

Cas. No. 3077385-98-6


Formula $C_{44}H_{45}F_4N_9O_5$ M.Wt 855.88

Solubility DMSO: 50 mg/mL (58.42 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

PROTAC STAT6 degrader-3 is a potent STAT6 PROTAC degrader with a DC_{50} of less than 1 nM. PROTAC STAT6 degrader-3 can be used for research on allergic/inflammatory diseases and cancer.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1].

In Vitro, PROTAC STAT6 degrader-3 is a potent STAT6 PROTAC degrader with a IC_{50} of less than 1 nM^[1].

References:

[1]. Bin Yang, et al. STAT6 DEGRADERS AND USES THEREOF CROSS REFERENCE TO RELATED APPLICATIONS. WO2025049820A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA