

Product Data Sheet

Product Name: Suc-Ala-Ala-Pro-Ala-pNA

Cat. No.: GA23527

Chemical Properties

Cas. No. 72682-69-0

Formula $C_{24}H_{32}N_6O_9$ M.Wt 548.55

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Suc-AAPA-pNA, an excellent substrate for porcine pancreatic elastase. Substrate for α -lytic protease.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA