

Product Data Sheet

Product Name: Tropisetron

Cat. No.: GC16284

Chemical Properties

Cas. No. 105826-92-4

Chemical Name [(1R,5S)-8-methyl-8-azabicyclo[3.2.1]octan-3-yl] 1H-indole-3-carboxylate

SMILES CN1C2CCC1CC(C2)OC(=O)C3=CNC4=CC=CC=C43


Formula C₁₇H₂₀N₂O₂ M.Wt 284.35

Solubility ≥ 28.4mg/mL in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Tropisetron (SDZ-ICS-930 free base) is a selective 5-HT₃ receptor antagonist and α₇-nicotinic receptor agonist with an IC₅₀ of 70.1 ± 0.9 nM for 5-HT₃ receptor. IC₅₀ value: 70.1 ± 0.9 nM [1]Target: 5-HT₃ receptor in vitro: Tropisetron specifically inhibited both IL-2 gene transcription and IL-2 synthesis in stimulated T cells. tropisetron inhibited both the binding to DNA and the transcriptional

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

activity of NFAT and AP-1. We also observed that tropisetron is a potent inhibitor of PMA plus ionomycin-induced NF-(kappa)B activation but in contrast TNF(alpha)-mediated NF-(kappa)B activation was not affected by this antagonist [2]. Tropisetron prevents the phosphorylation and thus activation of the p38 MAPK, which is involved in post-transcriptional regulation of various cytokines [3].in vivo: Two different doses of tropisetron (5 and 10 mg/kg) or vehicle were administered intraperitoneally 30 min before pMCAO. Neurological deficit scores, mortality rate and infarct volume were determined 24 h after permanent focal cerebral ischemia [4].

References:

- [1]. Macor JE, et al. The 5-HT₃ antagonist tropisetron (ICS 205-930) is a potent and selective alpha₇ nicotinic receptor partial agonist. *Bioorg Med Chem Lett*. 2001 Feb 12;11(3):319-21.
- [2]. Vega Lde L, et al. The 5-HT₃ receptor antagonist tropisetron inhibits T cell activation by targeting the calcineurin pathway. *Biochem Pharmacol*. 2005 Aug 1;70(3):369-80.
- [3]. Stratz C, et al. The anti-inflammatory effects of the 5-HT receptor antagonist tropisetron are mediated by the inhibition of p38 MAPK activation in primary human monocytes. *Int Immunopharmacol*. 2012 Aug;13(4):398-402.
- [4]. Candelario-Jalil E, et al. Detrimental effects of tropisetron on permanent ischemic stroke in the rat. *BMC Neurosci*. 2008 Feb 6;9:19.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA