

Product Data Sheet

Product Name: NVS-PI3-4

Cat. No.: GC62141

Chemical Properties

Cas. No. 941580-60-5


Formula $C_{20}H_{26}N_4O_3S$ M.Wt 402.51

Solubility DMSO : 250 mg/mL (621.10 mM; Storage Store at -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

NVS-PI3-4 is a specific PI3Ky inhibitor. NVS-PI3-4 can be used for the research of allergies,

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

inflammatory and cancer diseases[1][2].

NVS-PI3-4 shows an exquisite cellular selectivity for PI3K γ . NVS-PI3-4 reduces IgE/antigen-mediated phosphorylation of PKB/Akt in p110 δ DA. NVS-PI3-4 (5 μ M; 30 minutes; BMMCs) is not accumulate in specifically in mast cells[2].

[1]. Bruce I, et al. Development of isoform selective PI3-kinase inhibitors as pharmacological tools for elucidating the PI3K pathway. Bioorg Med Chem Lett. 2012;22(17):5445-5450.

[2]. Collmann E, et al. Transient targeting of phosphoinositide 3-kinase acts as a roadblock in mast cells' route to allergy. J Allergy Clin Immunol. 2013;132(4):959-968.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA