

## Product Data Sheet

Product Name: CAP1-6D

Cat. No.: GC74421

### Chemical Properties

Cas. No.


Formula  $C_{43}H_{68}N_{10}O_{15}$  M.Wt 965.06

Solubility Storage Sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



### Background

CAP1-6D is an agonist peptide for cytotoxic T cell (CTL), that enhances the immunogenicity of the CAPI peptide, and stimulates T cell responses.

#### References:

[1]. Zaremba S, et al., Identification of an enhancer agonist cytotoxic T lymphocyte peptide from human carcinoembryonic antigen. Cancer Res. 1997 Oct 15;57(20):4570-7.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA