

Product Data Sheet

Product Name: (Z)-4-(2-Cyano-2-phenylvinyl)benzonitrile

Cat. No.: GF48050

Chemical Properties

Cas. No. 61469-71-4

Formula C₁₆H₁₀N₂

M.Wt 230.26

Solubility

Storage 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA