

Product Data Sheet

Product Name: Copper chloride dihydrate AR, 99.0%

Cat. No.: GB57819

Chemical Properties

Cas. No. 10125-13-0

Formula $\text{Cl}_2\text{CuH}_4\text{O}_2$

M.Wt 170.48

Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA