

Product Data Sheet

Product Name: 2-Methyl-4-(piperidin-1-yl)phenylboronic acid

Cat. No.: GF34756

Chemical Properties

Cas. No. 2159136-20-4

Formula C₁₂H₁₈BN₂O₂

M.Wt 219.08782

Solubility

Storage -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA