

## Product Data Sheet

Product Name: 2-[3, 5-bis(trifluoromethyl)phenyl]-2-methylpropanoic acid

Cat. No.: GF38267

### Chemical Properties

Cas. No. 289686-70-0

Formula C<sub>12</sub>H<sub>10</sub>F<sub>6</sub>O<sub>2</sub>

M.Wt 300.19

Solubility

Storage Store at RT, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA