

Product Data Sheet

Product Name: LL-37 scrambled peptide

Cat. No.: GC36467

Chemical Properties

Cas. No.

Formula C₂₀₅H₃₄₀N₆₀O₅₃ M.Wt 4493.3

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

GLKLRFEFSKIKGEFLKTPEVRFDRDIKLDNRISVQR

Background

LL-37 scrambled peptide is a scrambled version of cathelicidin anti-microbial peptide LL-37. Scrambled peptide is random permutation of the original peptide that can be used as negative control or screening tool in research.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA