

Product Data Sheet

Product Name: GAPDH siRNA Positive Control (Mouse/Rat)

Cat. No.: GC81660

Chemical Properties

Cas. No.

Formula

M.Wt

15206.29

Solubility

Storage

Store at -20°C, away from moisture

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition

Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

GAPDH siRNA Positive Control (Mouse/Rat), a siRNA product, is a mouse/rat GAPDH siRNA positive control.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA