

Product Data Sheet

Product Name: L-R4W2 acetate

Cat. No.: GC79531

Chemical Properties

Cas. No.


Formula $C_{46}H_{71}N_{21}O_6 \cdot xC_2H_4O_2$ M.Wt 1014.19

Solubility Storage Sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Background

L-R4W2 acetate is a potent antagonist of vanilloid receptor 1 (VR1, TRPV1), with an IC₅₀ of 0.1 μ M.

L-R4W2 acetate may act as a potent analgesic [1] [2].

References:

[1]. Planells-Cases R, et, al. Arginine-rich peptides are blockers of VR-1 channels with analgesic activity. FEBS Lett. 2000 Sep 15; 481(2): 131-6.

[2]. Himmel HM, et, al. The arginine-rich hexapeptide R4W2 is a stereoselective antagonist at the vanilloid receptor 1: a Ca²⁺ imaging study in adult rat dorsal root ganglion neurons. J Pharmacol Exp Ther. 2002 Jun; 301(3): 981-6.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA