

Product Data Sheet

Product Name: Tetraammineplatinum(II) chloride monohydrate

Cat. No.: GD23098

Chemical Properties

Cas. No. 13933-33-0

Formula $\text{Cl}_2\text{H}_{14}\text{N}_4\text{OPt}$

M.Wt 352.121


Solubility

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA