

Product Data Sheet

Product Name: Boc-Ser(PO₃Bzl₂)-OH

Cat. No.: GA10513

Chemical Properties

Cas. No. 90013-45-9

Chemical Name (2S)-3-bis(phenylmethoxy)phosphoryloxy-2-[(2-methylpropan-2-yl)oxycarbonylamino]propanoic acid

SMILES CC(C)(C)OC(=O)NC(COP(=O)(OCC1=CC=CC=C1)OCC2=CC=CC=C2)C(=O)O

Formula C₂₂H₂₈NO₈P

M.Wt 465.3


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-Ser(PO₃Bzl₂)-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA