

## Product Data Sheet

---

Product Name: Lonigutamab

Cat. No.: GC26473

### Chemical Properties

Cas. No.

Formula

M.Wt

Solubility

Storage      Store at -80°C

General tips      For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition      Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

# Lonigutamab

### Background

Lonigutamab (hz208F2-4) is a humanized anti-IGF-1R monoclonal antibody that can be used for the

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

## Product Data Sheet

---

synthesis of antibody–agent conjugates (ADC)<sup>[1]</sup>.

### References:

[1]. Akla B, et al. Efficacy of the Antibody-Drug Conjugate W0101 in Preclinical Models of IGF-1 Receptor Overexpressing Solid Tumors. Mol Cancer Ther. 2020 Jan;19(1):168-177.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: [tech@glpbio.com](mailto:tech@glpbio.com)

Address: 10292 Central Ave. #205, Montclair, CA, USA