

Product Data Sheet

Product Name: L-Lysine thioctate

Cat. No.: GC34648

Chemical Properties

Cas. No. 20902-53-8

SMILES N[C@@H](CCCCN)C(O)=O.O=C(O)CCCCC1SSCC1


Formula $C_{14}H_{28}N_2O_4S_2$ M.Wt 352.51

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

L-Lysine thioctate is a substrate of lipoamidase[1].

[1]. Motafakkerazad R, et al. Simple HPLC evaluation of lipoamidase activity in tissue using a newly synthesized fluorescent substrate, dansyl- α -lipoyllysine. J Nutr Sci Vitaminol (Tokyo). 2011;57(5):377-82. Erratum in: J Nutr Sci Vitaminol (Tokyo). 2012;58(1):68.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA