

Product Data Sheet

Product Name: (Rac)-MEM 1003

Cat. No.: GC61873

Chemical Properties

Cas. No. 165187-25-7

SMILES O=C(OCCOC)C1=C(NC(C)=C(C1C2=CC=CC(C#N)=C2Cl)C(OC(C)C)=O)C


Formula $C_{22}H_{25}ClN_2O_5$ M.Wt 432.9

Solubility DMSO : 125 mg/mL (288.75 mM; ultrasonic and warming and heat to 80°C) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

(Rac)-MEM 1003 is the racemate of MEM 1003. MEM 1003, a dihydropyridine compound, is a potent L-type Ca^{2+} channel antagonist and has the potential for Alzheimer's disease research[1].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Gregory M Rose, et al. Efficacy of MEM 1003, a novel calcium channel blocker, in delay and trace eyeblink conditioning in older rabbits. Neurobiol Aging. 2007 May;28(5):766-73.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA