

Product Data Sheet

Product Name: PF-06737007

Cat. No.: GC70342

Chemical Properties

Cas. No. 1863905-38-7


Formula $C_{25}H_{28}F_4N_2O_6$ M.Wt 528.49

Solubility DMSO : ≥ 100 mg/mL (189.22 mM) Storage -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

PF-06737007 is a potent pan-Trk inhibitor in cell-based assays with IC_{50} s of 7.7 nM, 15 nM and 3.9 nM for TrkA, TrkB and TrkC, respectively. Anti-hyperalgesic effect.

PF-06737007 (Compound 13b) exhibits superb Trk selectivity with >95% inhibition of TrkA with no other kinase being inhibited by >40% at $1\text{ }\mu\text{M}$ ^[1].

References:

[1]. Bagal SK, et al. Discovery of Potent, Selective, and Peripherally Restricted Pan-Trk Kinase Inhibitors for the Treatment of Pain. J Med Chem. 2018 Aug 9;61(15):6779-6800.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA