

Product Data Sheet

Product Name: Iron(III)phosphatehydrate

Cat. No.: GB58955

Chemical Properties

Cas. No. 13463-10-0

Formula $\text{FeH}_4\text{O}_6\text{P}$

M.Wt 186.85


Solubility

Storage 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA