

Product Data Sheet

Product Name: Thalidomide-O-PEG4-azide

Cat. No.: GC66642

Chemical Properties

Cas. No. 2380318-57-8


Formula $C_{23}H_{29}N_5O_9$ M.Wt 519.5

Solubility DMSO : 250 mg/mL (481.23 mM; Storage Store at -20°C, away from moisture
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Thalidomide-O-PEG4-azide is a synthesized E3 ligase ligand-linker conjugate that incorporates the

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Thalidomide based cereblon ligand and a linker used in PROTAC technology.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins^[1].

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA