

Product Data Sheet

Product Name: Fmoc-Serinol(tBu)

Cat. No.: GA10757

Chemical Properties

Cas. No. 198561-87-4

Chemical Name (R,E)-(9H-fluoren-9-yl)methyl hydrogen (1-(tert-butoxy)-3-hydroxypropan-2-yl)carbonimide

SMILES CC(C)(OC[C@H](/N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O)([H])CO)C

Formula $C_{22}H_{27}NO_4$

M.Wt 369.4


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Serinol(tBu)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA