

Product Data Sheet

Product Name: N-Benzylaniline

Cat. No.: GD18421

Chemical Properties

Cas. No. 103-32-2

Formula

M.Wt


Solubility DMSO : 100 mg/mL (545.70 mM; Need ultrasonic)

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

N-Benzylaniline (N-Phenylbenzylamine) is an N-alkylated derivative of aniline. N-benzyl aniline also is a major metabolite of the antihistamine antazoline and other N-substituted benzyl anilines. N-Benzylaniline can be used for the research of various biochemical studies.

[1]. R Betz, et al. N-Benzylaniline. Acta Cryst. (2011). E67, o1195

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA