

Product Data Sheet

Product Name: GSK962
Cat. No.: GC38705

Chemical Properties

Cas. No. 2049872-86-6

SMILES CC(C)(C)C(N1N=CC[C@@H]1C2=CC=CC=C2)=O


Formula $C_{14}H_{18}N_2O$ M.Wt 230.31

Solubility DMSO: 75 mg/mL (325.65 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

GSK962 is an inactive enantiomer of GSK963 and can be used to confirm on-target effects[1].

[1]. Berger SB, et al. Characterization of GSK'963: a structurally distinct, potent and selective inhibitor of RIP1 kinase. Cell Death Discov. 2015 Jul 27;1:15009.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA