

Product Data Sheet

Product Name: 3,8-Bis[3,5-bis(trifluoromethyl)phenyl]-1,10-phenanthroline

Cat. No.: GD05030

Chemical Properties

Cas. No. 1228032-35-6

Formula C₂₈H₁₂F₁₂N₂

M.Wt 604.389


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA