

## Product Data Sheet

Product Name: ORIC-944 TFA

Cat. No.: GC73785

### Chemical Properties

Cas. No.

Formula  $C_{28}H_{26}F_4N_6O_3$

M.Wt 570.54


Solubility

Storage 4°C, sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



### Background

ORIC-944 TFA is the TFA salt form of ORIC. ORIC-944 TFA is a selective, orally active, allosteric inhibitor targeting the EED subunit of polycomb repressive complex 2 (PRC2). ORIC-944 TFA is synergistic with androgen receptor pathway inhibitors (ARPIs) for the study of metastatic prostate cancer.

References:

[1]. Daemen A, et al. ORIC-944, a potent and selective allosteric PRC2 inhibitor with best-in-class

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

## Product Data Sheet

---

properties, demonstrates combination synergy with AR pathway inhibitors in prostate cancer models[J]. Cancer Research, 2024, 84(6\_Supplement): 6586-6586. [2]. ORIC Pharmaceuticals Provides Initial Phase 1b Data for ORIC-944, Operational Highlights for 2023, and Anticipated Upcoming Milestones

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: [tech@glpbio.com](mailto:tech@glpbio.com)

Address: 10292 Central Ave. #205, Montclair, CA, USA